

#### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

13 Feb., 2003

# **Endothelium Preserving Microwave Treatment for Atherosclerosis**

Patrick Fink, Ph.D.

Arndt, G. D., Ph.D.

Ngo, Phong

NASA Johnson Space Center

**Avionic Systems Division** 



## **Background**

#### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Our primary focus is on antennas and electromagnetic systems
- We recently developed and licensed a microwave catheter antenna for cardiac arrhythmias
- We have developed catheter antennas and a system intended for treatment of atherosclerosis



## Major Technical Challenges

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Efficient transfer of energy through a small catheter
- High frequencies (short wavelength) due to restricted artery region
- Small, efficient antenna with directionality
- Verification of thermal profile
  - numerical simulation
  - relies on accurate knowledge of material properties
- Clinical testing




# **Technology Description**

#### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

13 Feb., 2003



Microwave catheter antenna for heating atherosclerotic lesions to reduce constriction in the artery

Antenna directionality focuses energy to lesions



## R&D Status

#### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Current level of development
  - Prototyped 2 directional catheter antennas specifically for atherosclerotic treatment
    - Open-ended waveguide with bevel (patented)
    - New directional catheter antenna
      - in disclosure process
  - Preliminary computational electromagnetic modeling (CEM) of new directional catheter antenna
  - Preliminary simulations of thermal profiles using candidate frequencies, power levels, and heating times
  - Performed initial test of energy transfer to lipids




## New Catheter Antenna

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

13 Feb., 2003



## Initial CEM simulation

- Ez near feed region
- Simplified infinitesimal current filament source
  - Demonstrates front/back directionality



## **Development Hurdles**

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

13 Feb., 2003

 Comprehensive test program with partner in the medical community

 Verify preservation of intima while supplying sufficient heat to atherosclerotic region



# **Competing Technologies**

**AVIONIC SYSTEMS DIVISION** 

NASA JSC/EV41/P. Fink

- Balloon angioplasty
- stents
- Laser ablation



## Comparison with Existing Technologies

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Disadvantages of current methods
  - ◆ Balloon angioplasty
    - ✓ Restenosis
  - ◆ stents
    - ✓ Restenosis in older, non-coated stents
  - Laser ablation
    - Possible destruction of artery



## Comparison with Existing Technologies

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Advantages of proposed method
  - We anticipate this method will be characterized by:
    - Prevention of restenosis
    - Inexpensive catheterization procedure on an outpatient basis
  - We are also extending this procedure to treat older, non-coated stents



## **Future R&D**

### **AVIONIC SYSTEMS DIVISION**

NASA JSC/EV41/P. Fink

- Comprehensive test program with partner in the medical community
- Future enhancement/extension:
  - possible extension to stents
  - ◆ Initiated test with stent in phantom material
- Remaining milestones
  - ◆ Test in extracted tissue
  - High fidelity computational model



## **IP Status**

**AVIONIC SYSTEMS DIVISION** 

NASA JSC/EV41/P. Fink

13 Feb., 2003

5 U.S. patents issued

- NASA Partnership Options
  - Space Act Agreement
    - Cooperative effort to mature the technology
  - Eventual licensing of technology